

Part III Patient Encounters and Performance Skills (PEPS®) Examination

SITE INFORMATION & EQUIPMENT LIST

For candidates testing beginning August 2025

EDISON BUILDING, 2nd FLOOR
7910 MICROSOFT WAY
CHARLOTTE, NC 28273
800.969.3926, EXTENSION 2

EXAM DATES:
August 2025 - May 2026

The information listed in this Part III PEPS Examination Site Information and Equipment List is for informational purposes only. There may be additional equipment or supplies on-site during the exam that may not be reflected in this manual. **Additionally, consumable supplies (e.g., alcohol wipes, tissues, cotton-tipped applicators) may change without warning.**

This information is not intended to prepare candidates for the exam, but to provide a general overview of equipment and supplies used during the exam. Candidates are reminded to thoroughly review the evaluation forms and Candidate Guide to adequately prepare for the Part III PEPS Examination.

Part III PEPS EQUIPMENT QUICK REFERENCE:

All Stations:

- Candidate device will be a laptop; candidate will bring into each room

Patient Encounter Stations:

- Note board (Tagau LCD Writing Tablet 15 inch)

Anterior Segment Station:

- Reliance 6200-L Chair
- 7900-IC Floor Stand
- Haag-Streit 900 BI Slit Lamp Video System
- Goldmann applanation tonometer with tonometer prism
- Volk 4-Mirror Glass Gonio Fundus Lens with Flange


Posterior Segment Station:

- Reliance 6200-L Chair
- 7900-IC Floor Stand
- Haag-Streit 900 BI Slit Lamp Video System
- Keeler BIO Vantage Plus Slimline Wireless Plus Indirect System with +2.00 D Oculars
- Volk 20D Condensing Lens
- Volk Digital Series High Mag Lens
- Volk Digital Series Wide Field Lens
- Volk 78D & 90D Lenses

CANDIDATE DEVICE

Candidates must log into the device when directed. Candidates will have the device during orientation and throughout the exam. Please review the tutorial on our website for specific instructions on using the platform during the exam.

For illustration purposes only, may not be exact representation




During the exam, the candidate will be able to use their device to:

- View station overview information
- View patient encounter clinical data
- Order and review ancillary tests
- Complete assessment and plans
- Submit incidents
- Monitor break, observation/overview, and exam cycle time

PATIENT ENCOUNTER (PE) STATIONS ROOM OVERVIEW (WITH CANDIDATE DEVICE SHOWN)

For illustration purposes only, may not be exact representation



PATIENT ENCOUNTER (PE) STATIONS

Each PE Station includes a charging cable and a mouse for optional use.



If a candidate wishes to use the mouse, candidates will need to connect the device. When connecting, ensure the cable is fully pushed into the port, as displayed in the image to the left. If you have difficulty connecting, inform the SP or the hallway monitor.

NOTE BOARD



Candidates may take notes in the PE rooms during observation time and throughout the exam cycle using a note board.

Note boards must remain in the exam room.

The note board uses a stylus to write. Slide out the stylus to use.

Lock and Unlock

The note board will be in the locked position when the candidate enters the room. The locked position will prevent the candidate from accidentally clearing their notes. If a candidate wishes to clear their notes, they must unlock it. The unlocking toggle switch is on the back of the device.



Clear

To clear out all data on the board, unlock the board and press the button on the right.

Once cleared, the previous information can not be restored.

SKILLS STATIONS

LIGHTING

Please check back for more information on room lighting

EXAM CHAIR

Exam Chair Overview Reliance 6200-L Chair



Exam Chair Overview


Reliance 6200-L Chair

To recline the exam chair:

Pull DOWN on the lever located on the back of the chair

Lower the chair and position appropriately

Release the lever while lowered to lock the chair in place



To rotate the chair, push the silver handle towards the floor. When the chair is in position, pull handle to the original position to lock in place.

SKILLS STATIONS

SLIT LAMP OVERVIEW


Slit Lamp

Haag-Streit 900 BI Slit
Lamp Video System


CAMERA IS ON THE LEFT OCULAR

For optimal views,
candidates must use
focusing rod to focus
the oculars

Bright White:
Top pin fixed all the way
to the left,
cannot be moved




Magnification:
10X or 16X
(left) (right)




Power:
blinking green = **standby**
solid green = **on**

Non-functioning buttons



Rheostat





Blue Filter

Circular dial turned all the way to the **left** so "B" is displayed



Yellow Filter

If the lever is at , the filter is IN

If the lever is at , the filter is OUT

(the filter cannot be in between these positions)



ANTERIOR SEGMENT SKILLS STATION ROOM OVERVIEW

For illustration purposes only, may not be exact representation



ANTERIOR SEGMENT SKILLS STATION


ADDITIONAL EQUIPMENT AND SUPPLIES

Example of anterior segment station supplies set up on the counter ready for use



Tray

- Tray will contain ready to use, already-disinfected, gonio lens and tonometer probe
- Please keep tray closed until the start of the exam
- After using the instruments, please place on a clean paper towel



Gonio Lens

Large Volk 4-Mirror Glass
Gonio Fundus Lens with Flange



ANTERIOR SEGMENT SKILLS STATION

ADDITIONAL SUPPLIES

Repeat Policy

Posted on the **skills station** exam room door



Miscellaneous Supplies


Alcohol Wipes



Kim Wipes



Tissues




Saline Solution



Cotton Tip Applicators



Elbow Rest



POSTERIOR SEGMENT SKILLS STATION ROOM OVERVIEW

For illustration purposes only, may not be exact representation



POSTERIOR SEGMENT SKILLS STATION

BIO

Keeler BIO Vantage Plus Slimline Wireless Plus Indirect System

- The BIO headset will be hanging on the stand
- Without an appropriate target and working distance (18-25 inches) the target will appear decentered
- The BIO has been equipped with our camera system to ensure captured views




FIXED

- aperture
- filter
- brightness
- circle beam height
- +2.00 oculars

ADJUSTABLE

- fit behind head
- fit on top of head
- PD
- Pantoscopic tilt: angle of ocular to eyes



POSTERIOR SEGMENT SKILLS STATION

ADDITIONAL EQUIPMENT AND SUPPLIES

Lens Overview

The following lenses are available for candidate use during posterior evaluation:

- Volk 20D Lens
- Volk Digital Series High Mag Lens
- Volk Digital Series Wide Field Lens
- Volk 78D & 90D Lenses



Miscellaneous Supplies


Alcohol Wipes



Kim Wipes




Tissues



Repeat Policy

Posted on the **skills station** exam room door



Elbow Rest



NBEO would like to extend its collective best wishes for your success.